

# Final Report

## **Pneumonia Detection Challenge**

By,  
Abdul, Abhishek, Tanu, Manish, Ashutosh

# Agenda

- Summary of problem statement, data and findings
- Overview of the final process
- Step-by-step walk through of the solution
- Model evaluation
- Comparison to benchmark
- Visualization(s)
- Implications
- Limitations
- Closing Reflections

# Problem Statement

- Build a deep learning model to detect and classify images of potential pneumonia cases so that doctors can prioritize and expediate the review.
- The model should predict whether pneumonia exists in a given image.
- The application should also predict by predicting bounding boxes around areas of the lung.
  - Samples without bounding boxes are negative and contain no definitive evidence of pneumonia.
  - Samples with bounding boxes indicate evidence of pneumonia

# DataSet Used

We are using the datasets provided by RSNA on Kaggle, we have downloaded a copy of these files

- stage\_2\_train.csv - the training set. Contains patientIds and bounding box / target information.
- stage\_2\_detailed\_class\_info.csv - provides detailed information about the type of positive or negative class for each image.
- stage\_2\_train\_images – set of training images
- Stage\_2\_test\_images - set of test images

# Overview of the final process

## Exploratory Data Analysis-

Preliminary EDA insights -class\_info dataset


- 30227 rows and 2 columns
- patientId and class

- Missing values –nil
- 26684 unique patients
- Class distribution



```
[ ] classes
```

➤ No Lung Opacity / Not Normal	11821
Lung Opacity	9555
Normal	8851
Name: class, dtype: int64	

	patientId	class
0	0004cfab-14fd-4e49-80ba-63a80b6bddd6	No Lung Opacity / Not Normal
1	00313ee0-9eaa-42f4-b0ab-c148ed3241cd	No Lung Opacity / Not Normal
2	00322d4d-1c29-4943-afc9-b6754be640eb	No Lung Opacity / Not Normal
3	003d8fa0-6bf1-40ed-b54c-ac657f8495c5	Normal
4	00436515-870c-4b36-a041-de91049b9ab4	Lung Opacity



### Training - Architecture



### Production - Architecture

## Preliminary EDA insights -train labels dataset

- 30227 rows and 6 columns

- patientId and x, y coordinates, width, height of the bounding boxes and target indicating Pneumonic case or not

- Missing values –nil

- 26684 unique patients

	patientId	x	y	width	height	Target
0	0004cfab-14fd-4e49-80ba-63a80b6bddd6	NaN	NaN	NaN	NaN	0
1	00313ee0-9eaa-42f4-b0ab-c148ed3241cd	NaN	NaN	NaN	NaN	0
2	00322d4d-1c29-4943-afc9-b6754be640eb	NaN	NaN	NaN	NaN	0
3	003d8fa0-6bf1-40ed-b54c-ac657f8495c5	NaN	NaN	NaN	NaN	0
4	00436515-870c-4b36-a041-de91049b9ab4	264.0	152.0	213.0	379.0	1

## Preliminary EDA insights –training images

- Training has 26684 images
- Images are dcm format
- Every image has the following details
  - Patient sex;
  - Patient age;
  - Modality;
  - Body part examined;
  - View position;
  - Rows & Columns;
  - Pixel Spacing.

(0008, 0005) Specific Character Set	CS: 'ISO_IR 100'
(0008, 0016) SOP Class UID	UI: Secondary Capture Image
(0008, 0018) SOP Instance UID	UI: 1.2.276.0.7230010.3.1.4
(0008, 0020) Study Date	DA: '19010101'
(0008, 0030) Study Time	TM: '000000.00'
(0008, 0050) Accession Number	SH: ''
(0008, 0060) Modality	CS: 'CR'
(0008, 0064) Conversion Type	CS: 'WSD'
(0008, 0090) Referring Physician's Name	PN: ''
(0008, 103e) Series Description	LO: 'view: PA'
(0010, 0010) Patient's Name	PN: '0a2c130c-c536-4651-836
(0010, 0020) Patient ID	LO: '0a2c130c-c536-4651-836
(0010, 0030) Patient's Birth Date	DA: ''
(0010, 0040) Patient's Sex	CS: 'M'
(0010, 1010) Patient's Age	AS: '30'
(0018, 0015) Body Part Examined	CS: 'CHEST'
(0018, 5101) View Position	CS: 'PA'
(0020, 000d) Study Instance UID	UI: 1.2.276.0.7230010.3.1.2
(0020, 000e) Series Instance UID	UI: 1.2.276.0.7230010.3.1.3
(0020, 0010) Study ID	SH: ''
(0020, 0011) Series Number	IS: "1"
(0020, 0013) Instance Number	IS: "1"
(0020, 0020) Patient Orientation	CS: ''
(0028, 0002) Samples per Pixel	US: 1
(0028, 0004) Photometric Interpretation	CS: 'MONOCHROME2'
(0028, 0010) Rows	US: 1024
(0028, 0011) Columns	US: 1024
(0028, 0030) Pixel Spacing	DS: ['0.143000000000000002',
(0028, 0100) Bits Allocated	US: 8
(0028, 0101) Bits Stored	US: 8
(0028, 0102) High Bit	US: 7
(0028, 0103) Pixel Representation	US: 0
(0028, 2110) Lossy Image Compression	CS: '01'
(0028, 2114) Lossy Image Compression Method	CS: 'ISO_10918_1'
(7fe0, 0010) Pixel Data	OB: Array of 146200 bytes



## Algorithms Used-

YOLO/OpenCV DarkFlow –TensorFlow implementation of YOLO  
ChexNet to identify Pneumonia Patient

## Step-by-step walk through of the solution

- ❖ Build a deep learning model to detect and classify images of potential pneumonia cases so that doctors can prioritize and expedite their review.
- ❖ The model should predict whether pneumonia exists in a given image.
- ❖ The application should also predict by predicting bounding boxes around areas of the lung.
- ❖ Samples without bounding boxes are negative and contain no definitive evidence of pneumonia.
- ❖ Samples with bounding boxes indicate evidence of pneumonia.

# Model evaluation

- The evaluation metrics is based on the mean average precision at different intersection over union (IoU) thresholds. The IoU of a set of predicted bounding boxes and ground truth bounding boxes is calculated as:

$$\text{IoU}(A,B) = \frac{A \cap B}{A \cup B}$$

- The metric sweeps over a range of IoU thresholds, at each point calculating an average precision value. The threshold values range from 0.4 to 0.75 with a step size of 0.05: (0.4, 0.45, 0.5, 0.55, 0.6, 0.65, 0.7, 0.75). In other words, at a threshold of 0.5, a predicted object is considered a "hit" if its intersection over union with a ground truth object is greater than 0.5.

- At each threshold value  $t$ , a precision value is calculated based on the number of true positives (TP), false negatives (FN), and false positives (FP) resulting from comparing the predicted object to all ground truth objects:

$$\text{Precision}(t) = \frac{\text{TP}(t)}{\text{TP}(t) + \text{FP}(t)}$$

- A true positive is counted when a single predicted object matches a ground truth object with an IoU above the threshold. A false positive indicates a predicted object had no associated ground truth object. A false negative indicates a ground truth object had no associated predicted object.

**Important note:**if there are no ground truth objects at all for a given image, ANY number of predictions (false positives) will result in the image receiving a score of zero, and being included in the mean average precision.

The average precision of a single image is calculated as the mean of the above precision values at each IoUthreshold:

$$\frac{1}{|\text{thresholds}|} \sum_t \frac{\text{TP}(t)}{\text{TP}(t) + \text{FP}(t) + \text{FN}(t)}$$



- Bounding boxes will be evaluated in order of their confidence levels in the above process. This means that bounding boxes with higher confidence will be checked first for matches against solutions, which determines what boxes are considered true and false positives.



## Intersection over Union (IoU)

- Intersection over Union is a measure of the magnitude of overlap between two bounding boxes (or, in the more general case, two objects). It calculates the size of the overlap between two objects, divided by the total area of the two objects combined.

- The two boxes in the visualization overlap, but the area of the overlap is insubstantial compared with the area taken up by both objects together. IoU would be low -and would likely not count as a "hit" at higher IoU thresholds.

$$\textit{Intersection over Union (IoU)} = \frac{\textit{Area of Overlap}}{\textit{Area of Union}}$$

-  Prediction
-  Ground-truth



# Comparison to benchmark

- CheXNet accuracy is used as the standard for the benchmark purpose.
- Currently, accuracy achieved with this RSNA dataset using CheXNet is 76.8%
- We have created the model using YOLO and it is providing the accuracy as 78.2%


# Visualization

## Exploratory visualization

- Normal lungs are filled with air. In the x-rays below, you can see normal lungs. We can note that we mostly see white skeletal matter and black matter, which is primarily air.



When a person has pneumonia, however, air is replaced by fluid, bacteria, immune system cells, and other objects. In an xray, the opacities tend to have a greyish color, and a cloudy appearance, rather than being black or white. You can see some images with pneumonia related opacities below.






While we are theoretically detecting “lung opacities”, there are lung opacities that are not pneumonia related

In the data, some of these are confusingly labeled “Not Normal No Lung Opacity” These non pneumonia “Not Normal” detections end up being a primary source of frustration in building models

The images below show a few of the examples in the “Not Normal” class.



# Implications

We have build an algorithm to detect a visual signal for pneumonia in medical images to improve the efficiency and reach of diagnostic services.

**Here's the backstory and why solving this problem matters.**

Pneumonia accounts for over 15% of all deaths of children under 5 years old internationally. In 2017, 920,000 children under the age of 5 died from the disease. In the United States, pneumonia accounts for over 500,000 visits to emergency departments and over 50,000 deaths in 2017,, keeping the ailment on the list of top 10 causes of death in the country.

While common, accurately diagnosing pneumonia is a tall order. It requires review of a chest radiograph (CXR) by highly trained specialists and confirmation through clinical history, vital signs and laboratory exams. Pneumonia usually manifests as an area or areas of increased opacity on CXR. However, the diagnosis of pneumonia on CXR is complicated because of a number of other conditions in the lungs such as fluid overload (pulmonary edema), bleeding, volume loss (atelectasis or collapse), lung cancer, or post-radiation or surgical changes. Outside of the lungs, fluid in the pleural space (pleural effusion) also appears as increased opacity on CXR. When available, comparison of CXRs of the patient taken at different time points and correlation with clinical symptoms and history are helpful in making the diagnosis.

CXRs are the most commonly performed diagnostic imaging study. A number of factors such as positioning of the patient and depth of inspiration can alter the appearance of the CXR, complicating interpretation further. In addition, **clinicians are faced with reading high volumes of images every shift.**

# Limitations

1. Infrastructure - GPU
2. Quality labeled image
3. Access to radiology
  - Need to work with the hospitals to check the efficiency of the model

# Closing Reflections

- Shallow practical knowledge
- Process pipeline
- Domain understanding of health care side
- Team effort & collaboration
- Diversified profile of team mates.
- Good understanding of deep learning is required